

DIRECTIONS: Fill container with warm water to brim, secure cap and shake well. Repeat shaking-do not add additional water. This will produce a specific gravity of 1.18.

SHAKE WELL BEFORE USING

NDC 57319-360-09

OVUM FLOTATION DRY

Fecal Flotation Medium

Ovum Flotation Dry is a pre-measured amount of Zinc Sulfate to be reconstituted and used for flotation techniques of fecal analysis.

Protect from freezing.

Keep tightly closed when not in use.

TAKE TIME



OBSERVE LABEL DIRECTIONS

Net Contents: 1.32 Kg. (2.9 lbs.)

Lot No.

Exp. Date

Rev. 1-10

FOR ANIMAL USE ONLY

KEEP OUT OF REACH OF CHILDREN



Manufactured for:
Clipper Distributing Company, LLC
St. Joseph, MO 64507

Trademarks are property of
Clipper Distributing Company, LLC



3